

Medela Flex™
 technology
 The first pumping
 experience shaped

around you



PersonalFit Flex™ Breast Shields

PersonalFit™ Breast Shields are Medela's best-selling accessory*. With PersonalFit Flex, mom can easily adjust the way her breast shield fits to find the most comfortable pumping position, even as her breast shape changes. It's a whole new pumping experience, personalized for mom.

Benefits

- **Adaptable:** the new 105° angle and oval shape are designed to better fit the shape of mom's breast and help milk flow more easily.
- **Comfortable:** oval shape can be rotated 360° for comfortable, flexible positioning.
- **Effectiveness:** removes 11.8% more milk per minute.
- **Improved Efficiency/Faster Milk Ejection:** the new opening angle is designed to reduce pressure on the breast and optimize milk flow.

Features:

Oval shape can be rotated 360° for comfortable, flexible positioning



105° angle adapts to the breast



Soft, smooth, natural-feeling rim

Product Specifications:

Item #	Description	UPC Code	Qty/ Unit	Unit Dimensions (HxWxD) Inches	Qty Per Case Pack	Shipping Carton Dimensions (LxWxH) Inches
101036653	21mm PersonalFit Flex™ Breast Shield	020451366531	2	5.5 x 3.94 x 3.54	6	11.25 x 8.38 x 4.13
101036654	24mm PersonalFit Flex Breast Shield	020451366548	2	5.5 x 3.94 x 3.54	6	11.25 x 8.38 x 4.13
101036655	27mm PersonalFit Flex Breast Shield	020451366555	2	6.5 x 4.10 x 3.7	6	11.5 x 8.75 x 5.13
101036656	30mm PersonalFit Flex Breast Shield	020451366562	2	6.5 x 4.10 x 3.7	6	11.5 x 8.75 x 5.13

To order or for more information, contact your Medela Sales Consultant.