

# Medela Breast Pumps



## The #1 Recommended Breast Pump Brand\*



Pump in Style® Advanced Starter Set



Pump in Style® Advanced Backpack & On-the-go Tote



Pump In Style® Advanced The Metro Bag™



Freestyle Flex™ Breast Pump



Sonata® Breast Pump



Symphony® Breast Pump

SKU	101035077	101036453 / 101036449	101036445	101037980	58200	0240108 / 0240208
<b>Description Marketing or Packaging</b>	Medela's Pump in Style Advanced Starter Set delivers effective pumping in a modular solution. Offering trusted and research based performance with 2-Phase Expression® in a streamlined format customizable to mom's changing lifestyle.	Medela's award-winning Pump in Style Advanced breast pump designed for moms who need true portability when pumping on-the-go. Containing Medela's trusted research based performance of 2-Phase Expression® technology providing mom everything she needs to get started and stay organized while breast pumping.	The same great Pump in Style Advanced Breast Pump with a premium carry bag and extra accessories. This double electric breast pump, designed for portability, helps mom balance breast pumping with her goals & everyday life. It contains Medela's trusted research based performance of 2-Phase Expression® technology in an all-in-one breast pump system.	Freestyle Flex offers true mobility in a lightweight and compact double electric personal use breast pump that is small enough to fit in the palm of your hand or pocket. Designed for Life, Freestyle Flex offers Medela's original 2-Phase Expression® technology, the convenience of a USB-chargeable* 2-hour battery and touchscreen display with easy to use controls & timer. Freestyle Flex is perfect for the mom that prefers to multitask and move around while pumping. <small>*Only use the power adapter and cable that comes with the breast pump.</small>	Medela's Sonata Breast Pump offers our first connected solution to empower and support moms, breast milk feeding success through smart connectivity for confidence and convenience, innovative performance features designed to ensure optimal efficiency in every session and quiet operation for a comfortable experience.	Medela's Gold-Standard hospital-grade (multi-user) breast pump designed for breast milk initiation, building and maintenance.
<b>Technology</b>	2-Phase Expression	2-Phase Expression	2-Phase Expression	2-Phase Expression	2-Phase Expression	2-Phase Expression
<b>Daily Use</b>	●	●	●	●	●	●
<b>Variable Vacuum/Speed Settings</b>	● Combined & Fully Adjustable	● Combined & Fully Adjustable	● Combined & Fully Adjustable	● Combined & Fully Adjustable	● Combined & Fully Adjustable	● Combined & Fully Adjustable
<b>mmHg - Vacuum</b>	-50 to -240 mmHg, maximum of -295	-50 to -240 mmHg, maximum of -295	-50 to -240 mmHg, maximum of -295	-45 to -245 mmHg, maximum of -270	-40 to -250 mmHg	-30 to -270 mmHg
<b>Cycles Per Minute</b>	54-120 cpm	54-120 cpm	54-120 cpm	45 to 111 cpm	54-120 cpm	54-120 cpm
<b>Carry Bag</b>	Motor bag only (carry bag sold separately)	All-in-one portable black microfiber shoulder bag or backpack.	All-in-one portable black microfiber shoulder bag includes work surface and mesh bags for part organization.	All-in-one portable black microfiber shoulder bag	All-in-one portable black microfiber shoulder bag	Hard protective carry case
<b>Removable Cooler &amp; Ice Pack</b>	Sold Separately	●	●	●	●	Sold Separately
<b>Breast Milk Containers</b>	2 (5 oz)	4 (5 oz)	4 (5 oz)	4 (5 oz)	4 (5 oz)	2 (5 oz) included in the double pumping kit
<b>Battery</b>	8AA External Battery Pack (AA not included)	8AA External Battery Pack (AA not included)	8AA External Battery Pack (AA not included)	Rechargeable Lithium-Ion 2 hrs.	Rechargeable Lithium-Ion 1 hr.	Symphony Breast Pump with Rechargeable Battery 0240208 1 hr.
<b>Digital Display</b>	--	--	--	●	●	●
<b>Timer</b>	--	--	--	●	● Count-up and Count-down	--
<b>Hospital-Grade (multi-user)</b>	--	--	--	--	--	●
<b>Hospital Performance</b>	--	--	--	--	● Clinically Tested Hospital Performance Pumping Rhythms designed for maintaining milk supply	● Clinically Tested Hospital Performance Pumping Rhythms designed for initiation, building & maintaining milk supply
<b>Closed System with Overflow Protection</b>	● <sup>^</sup>	● <sup>^</sup>	● <sup>^</sup>	● <sup>#</sup>	● <sup>#</sup>	● <sup>#</sup>
<b>Bluetooth® Connectivity</b>	--	--	--	●	●	--

Content appearance may vary from picture.  
<sup>^</sup>Overflow Protection is located at the pump. <sup>#</sup>Overflow Protection is located at the kit.