






















Medela breast pump component compatibility chart

To ensure proper performance and avoid breast pump damage, use only the components designed to work with your breast pump.

Breast Pumps		Compatible Connector	Compatible Tubing	Part Number
	Pump In Style® WITH maxFLOW™		 Single Y- tubing : All three adaptor ends are narrow & round.	Connector 101041267 101038233 Tubing 101040485
	Freestyle Flex™		 Single Y- tubing : two narrow round ends and 1 larger triangle end.	Connector 101038223 101041267 Tubing 101038234
	Sonata®		 Single Y- tubing : two narrow round ends and 1 larger round end.	Connector 68054 Tubing 101038540
	Harmony®		NO TUBING	Connector 101034753 101034864
	Swing®		 One single 30" tube: with one round end and one triangle end.	Connector 101034753 Tubing 101034768
	Symphony®		 Two separate 42" tubes: with one round end and one triangle end.	Connector 87071 Tubing 101034773
	Pump In Style® Advanced		 Two separate tubes: with one round end and one soft tubing end.	Connector 101039190 87071 Tubing 101033078
	Freestyle™		 Single Y- tubing : All three adaptor ends are thick & round.	Connector 67061 Tubing 101033077