

Low Dose Dual Cap

DualCap meets ISO standard for syringe accuracy when used as indicated.


The Medela Dual Cap provides a safe, accurate, and convenient way for your pharmacy to continue the use of ISMP recommended slip-tip syringes, while allowing the entire hospital to benefit from transitioning to ENFit® enteral connectors. Medela's Dual Cap does not require twisting to apply the cap to the syringe. This feature results in no change to current filling practices, which may result in saved time and labor.

Benefits:

- Pharmacy can use the same oral slip-tip syringes and filling methods.
- 2-piece cap facilitates a clean connection.
- Removal of the first cap facilitates connection to ENFit feeding tubes at time of administration.
- Designed with Low Dose technology, which provides accurate dosing through the feeding tube.
- The Low Dose Dual Cap is only needed when administering medication through an ENFit feeding tube.



Just 3 easy steps....



Product Specifications:

Item #	Description	Sterile or Non-Sterile	Qty/Case (each)
ENFDUALCAPLD	ENFit to TwistLok® connector, female to female with cap	Sterile	500 (50 packs of 10) per case
ENFDUALCAPNSLD	ENFit to TwistLok connector, female to female with cap	Non-Sterile	500 caps per pouch

To order or for more information, contact your Medela Sales Consultant.

Content appearance may vary from picture.

Medela wordmark and logo, Every Drop Counts, and TwistLok are registered in the U.S. Patent and Trademark Office. ENFit is a federally registered trademark of GEDSA and is used with permission. Neonatal Solutions logo is a trademark of Medela.

1549104 B 0620 © 2020 Medela 4550174